
PCR EN TEMPO REAL

MELTING CURVE: CURVA DE DISOCIACIÓN

8/VII/2014

S1401014

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


Prememos **SETUP**.

PCR EN TEMPO REAL

ABI PRISM STEPONE: Experiment Properties

Pantalla de Experiment Properties

StepOne™ Software v2.1

File Edit Instrument Analysis Tools Help

New Experiment Open Save Close Send Experiment to Instrument Download Experiment from Instrument Print Print Report

Experiment: **Untitled** Type: **Quantitation - Standard Curve** Reagents: **TaqMan® Reagents**

Design your experiment...

1. Define

Experiment Properties

2. Set Up

3. Order (optional)

Materials List

1A. Define: Experiment Properties

Instructions: Enter an experiment name, select the instrument type, then select the type of experiment to design.

How do you want to identify this experiment? **Required**

Experiment Name: **Untitled**

Barcode (Optional):

User Name (Optional):

Comments (Optional):

Which instrument are you using to run the experiment?

StepOnePlus™ Instrument (96 Wells) **StepOne™ Instrument (48 Wells)**

Set up, run, and analyze an experiment using a 3-color, 48-well system.

What type of experiment do you want to design?

Quantitation Genotyping Presence / Absence

Design a gene quantitation experiment to determine the amount of target nucleic acid sequence in a sample.

Previous Finish Designing Experiment Next Cancel

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

PCR EN TEMPO REAL

ABI PRISM STEPONE: Methods & Materials

The screenshot displays the 'StepOne Software v2.3' interface. The main window is titled 'Design your experiment...' and shows the 'Methods & Materials' section for an experiment named 'Practica larvas'. The experiment type is 'Quantitation - Standard Curve' and the reagents are 'SYBR® Green Reagents'. The interface is divided into three main sections: 1. Define, 2. Set Up, and 3. Order (Optional). The '1. Define' section includes 'Experiment Properties', 'Methods & Materials', 'Targets', 'Standards', 'Samples', 'Run Method', and 'Reaction Setup'. The '2. Set Up' section includes 'Targets', 'Standards', 'Samples', 'Run Method', and 'Reaction Setup'. The '3. Order (Optional)' section includes 'Materials List'. The 'Methods & Materials' section contains the following questions and options:

- 1B. Define: Methods & Materials**
Instructions: Select the quantitation method, reagents, ramp speed, and type of template for the real-time PCR reactions.
- Which quantitation method are you using?**
Options: ☒ Standard Curve, ☐ Relative Standard Curve, ☐ Comparative Ct ($\Delta\Delta Ct$)
Note: With the standard curve method, you use standards to determine the absolute quantity of target sequence in a sample.
- Which reagents do you want to use to detect the target sequence?**
Options: ☐ TaqMan® Reagents, ☒ SYBR® Green Reagents
Note: The real-time PCR reactions contain two primers and SYBR® Green I dye. The primers are designed to amplify the target sequence. The SYBR Green dye generates increased fluorescence signal when it binds to double-stranded DNA. Select "Include Melt Curve" to perform melt curve analysis of the amplified sequence.
☒ Include Melt Curve
- Which ramp speed do you want to include in the instrument run?**
Options: ☒ Standard (~ 2 hours to complete a run), ☐ Fast (~ 40 minutes to complete a run)
Note: For optimal results using the standard ramp speed, Applied Biosystems recommends standard reagents for your real-time PCR reactions.
- What type of template do you want to use in the real-time PCR reactions?**
Options: ☐ cDNA (complementary DNA), ☐ RNA, ☒ gDNA (genomic DNA)
Note: You are adding purified gDNA to the real-time PCR reactions. You have already extracted the gDNA from tissue or sample.


At the bottom of the window, there are navigation buttons: 'Previous', 'Finish Designing Experiment', 'Next', and 'Cancel'. The status bar at the bottom indicates 'Local Instrument: connected'.

Curva Standard, curva estandard relativa, Comparativa de C_T
Reagents. Tipo de química:
SYBR Green, TaqMan

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

PCR EN TEMPO REAL

ABI PRISM STEPONE: SETUP STANDARD

The screenshot displays the StepOne Software v2.3 interface. The main window is titled '2B. Set Up: Standards'. It includes a sidebar on the left with navigation options: '1. Define', '2. Set Up', and '3. Order (Optional)'. The '2. Set Up' section is active, showing 'Targets', 'Standards', 'Samples', 'Run Method', and 'Reaction Setup'. The 'Standards' section is highlighted.

The 'Set Up Standards' panel on the left contains the following information:

- Instructions:** Enter the number of points and replicates for all standard curves in the reaction plate. For each standard curve, enter the starting quantity and select the serial factor.
- Set Up Standards:**
 - How many **points** do you need for each standard curve?
 - How many **replicates** do you need for each point?
 - For each standard curve, define the range of standard quantities by entering the **starting quantity** and **serial factor**.
 - Define the standard curve so that the range of standard quantities spans the expected range of quantities for the unknowns.
 - Table with columns: Target Name, Enter Starting Quantity, Select Serial Factor.
 - Row 1: 16s vpul, 1.0, 2*
- Standard Curve Preview:** A graph showing a standard curve for 16s vpul. The y-axis is labeled with values: 3.2E1, 1.6E1, 8E0, 4E0, 2E0, 1E0. The x-axis represents the serial factor.
- Well Count:** 6 - Unknown, 30 - Standard, 3 - Negative Control, 9 - Empty.

The 'View Plate Layout' panel on the right shows a 96-well plate layout. It includes options to 'Arrange Plate by: Columns' and 'Place Negative Controls in: Lower Right'. The plate layout is displayed as a grid with columns 1-8 and rows A-F. The layout shows the placement of standards, samples, and negative controls.

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

ABI PRISM STEPONE: REACTION SETUP

StepOne Software v2.3

File Edit Instrument Analysis Tools Help

New Experiment Open Save Send Experiment to Instrument Download Experiment from Instrument Print Report

Design your experiment... Experiment: **Practica larvas** Type: **Quantitation - Standard Curve** Reagents: **SYBR® Green Reagents**

1. Define

2. Set Up

3. Order (Optional)

Materials List

2F. Set Up: Reaction Setup

Instructions: For each target assay in the reaction plate, select the assay type (if using TaqMan reagents), then review the calculated volumes for preparing the standard dilution series, samples, and PCR reactions. If needed, edit the reaction volume, excess reaction volume, component concentrations, and/or stock concentrations. Click "Print Reaction Setup" to print instructions on how to prepare the PCR reactions.

Select Tab: **Reaction Mix Calculations** Sample Dilution Calculations

Reaction Volume Per Well: 13.0 µL Excess Reaction Volume: 10 %

Master Mix Concentration: 2.0 x

Forward Primer Starting Concentration: 20.0 pmol/µL Forward Primer Final Concentration: 300.0 nM

Reverse Primer Starting Concentration: 20.0 pmol/µL Reverse Primer Final Concentration: 300.0 nM

Component	Volume (µL) for 1 Reaction
Master Mix (2.0x)	6.50
Forward Primer	0.20
Reverse Primer	0.20
Template (1.0x or Standard)	1.30
H ₂ O	4.81
Total Volume	13.00

Standard Dilution Series for 16x vpu

Dilution Point	Source	Source Volume (µL)	Diluent Volume (µL)	Total Volume (µL)	Standard Concentration (ng/µL)
1 (C2)	Stock	3.87	11.86	15.73	24.6
2 (F8)	Dilution 1	7.86	7.86	15.73	12.3
3 (F8)	Dilution 2	7.86	7.86	15.73	6.2
4 (F8)	Dilution 3	7.86	7.86	15.73	3.1
5 (C2)	Dilution 4	7.86	7.86	15.73	1.6
6 (C1)	Dilution 5	7.86	7.86	15.73	0.8

Standard Concentration in Stock: 100.0 ng µL

Standard Concentration in Stock: 100.0 ng µL

Previous Finish Designing Experiment Next

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
New Experiment Open Save Close Send Experiment to Instrument Download Experiment from Instrument Report Print Report

Experiment Menu Experiment: **Practica larvas** Type: **Standard Curve** Reagents: **SYBR® Green Reagents** **START RUN**

Setup
Experiment Properties
Plate Setup
Run Method
Reaction Setup
Materials List
Run
Analysis

Reaction Setup

Instructions: For each target assay in the reaction plate, select the assay type (if using TapMan reagents), then review the calculated volumes for preparing the standard dilution series, samples, and real-time PCR reactions. If needed, edit the reaction volume, excess reaction volume, component concentrations, and/or stock concentrations. Click "Print Reaction Setup" to print instructions on how to prepare the PCR reactions.

Reaction Mix Calculations **Sample Dilution Calculations**

Select Tab: **Reaction Mix Calculations** Reaction Volume Per Well: 12 µL Excess Reaction Volume: 10 % **Print Reaction Setup**

Reactions for 16s vpu

Master Mix Concentration: 2.0 x
Forward Primer Starting Concentration: 20.0 pmol/µL Forward Primer Final Concentration: 300.0 nM
Reverse Primer Starting Concentration: 20.0 pmol/µL Reverse Primer Final Concentration: 300.0 nM

Component	Volume (µL) for 1 Reaction
Master Mix (2.00)	6.90
Forward Primer	0.18
Reverse Primer	0.18
Sample (1.00) or Standard	1.20
H ₂ O	4.44
Total Volume	12.00

Standard Dilution Series for 16s vpu

Dilution Point	Source	Source Volume (µL)	Diluent Volume (µL)	Total Volume (µL)	Standard Concentration (ng/µL)
1 (32)	Stock	3.87	18.65	14.52	26.7
2 (16)	Dilution 1	7.26	7.26	14.52	13.3
3 (8)	Dilution 2	7.26	7.26	14.52	6.7
4 (4)	Dilution 3	7.26	7.26	14.52	3.3
5 (2)	Dilution 4	7.26	7.26	14.52	1.7
6 (1)	Dilution 5	7.26	7.26	14.52	0.8

Starting run
Starting the instrument run...

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New Experiment Open Save Close Send Experiment to Instrument... Download Experiment from Instrument... Export Print Report

Experiment Menu Experiment: **Practica larvas** Type: **Standard Curve** Reagents: **SYBR® Green Reagents** **START RUN**

Setup

- Experiment Properties
- Plate Setup
- Run Method
- Reaction Setup**
- Materials List

Run

Analysis

Reaction Setup

Instructions: For each target assay in the reaction plate, select the assay type (if using TaqMan reagents), then review the calculated volumes for preparing the standard dilution series, samples, and real-time PCR reactions. If needed, edit the reaction volume, excess reaction volume, component concentrations, and/or stock concentrations. Click "Print Reaction Setup" to print instructions on how to prepare the PCR reactions.

Reaction Mix Calculations **Sample Dilution Calculations**

Select Target: **16s vpul**

Reaction Volume Per Well: **12** μ L Excess Reaction Volume: **10** %

Reactions for 16s vpul

Master Mix Concentration: **2.0** X

Forward Primer Starting Concentration: **20.0** pmol/ μ L Forward Primer Final Concentration: **300.0** nM

Reverse Primer Starting Concentration: **20.0** pmol/ μ L Reverse Primer Final Concentration: **300.0** nM

Component	Volume (μ L) for 1 Reaction
Master Mix (2.0X)	6.00
Forward Primer	0.18
Reverse Primer	0.18
Sample (10X) or Standard	1.20
H ₂ O	4.44
Total Volume	12.00

Starting run

Starting the instrument run...

Standard Dilution Series for 16s vpul

Standard Concentration in Stock: **100.0** ng per μ L


Dilution Point	Source	Source Volume (μ L)	Diluent Volume (μ L)	Total Volume (μ L)	Standard Concentration (ng/ μ L)
1 (32)	Stock	3.87	10.65	14.52	26.7
2 (16)	Dilution 1	7.26	7.26	14.52	13.3
3 (8)	Dilution 2	7.26	7.26	14.52	6.7
4 (4)	Dilution 3	7.26	7.26	14.52	3.3
5 (2)	Dilution 4	7.26	7.26	14.52	1.7
6 (1)	Dilution 5	7.26	7.26	14.52	0.8

Local Instrument: connected.

1:45 PM

PCR EN TEMPO REAL

RUN

PCR EN TEMPO REAL

RUN